
Product Datasheet

KLHL15 Antibody (orb1244234)

Description

KLHL15 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody used in WB on Human 293T at 0.2...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human KLHL15.

Target

KLHL15

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store KLHL15 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

KLHL15 antibody can be used for detection of KLHL15 by ELISA at 1:312500. KLHL15 antibody can be used for detection of KLHL15 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

66 kDa

Uniprot ID
[Q96M94](#)
NCBI
[NP_085127](#)
Dilution Range

KLHL15 antibody can be used for detection of KLHL15 by ELISA at 1:312500. KLHL15 antibody can be used for detection of KLHL15 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.