



Product Datasheet

ERLIN2 Antibody (orb1244237)

Description	ERLIN2 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding to a region of human ERLIN2.	Antibody used in WB on Human 293T at 0.2...
Target	ERLIN2	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store ERLIN2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	Antibody used in WB on Human HeLa, AT3 a...
Note	For research use only	
Application notes	ERLIN2 antibody can be used for detection of ERLIN2 by ELISA at 1:62500. ERLIN2 antibody can be used for detection of ERLIN2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	38 kDa	
Uniprot ID	O94905	Antibody used in WB on Hum. Adult Placen...
NCBI	NP_009106	
Dilution Range	ERLIN2 antibody can be used for detection of ERLIN2 by ELISA at 1:62500. ERLIN2 antibody can be used for detection of ERLIN2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	