

Product Datasheet


SPATA7 Antibody (orb1244303)

Description


SPATA7 Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SPATA7.
Target	SPATA7
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store SPATA7 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	SPATA7 antibody can be used for detection of SPATA7 by ELISA at 1:1562500. SPATA7 antibody can be used for detection of SPATA7 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	68 kDa
Uniprot ID	Q9P0W8
NCBI	NP_060888
Dilution Range	SPATA7 antibody can be used for detection of SPATA7 by ELISA at 1:1562500. SPATA7 antibody can be used for detection of SPATA7 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.

Antibody used in WB on Human 293T at 0.2...



Antibody used in WB on Hum. Adult Placen...



Antibody used in WB on Hum. Fetal Heart ...