
Product Datasheet

ALKBH3 Antibody (orb1244386)

Description

ALKBH3 Antibody

Species/Host

Rabbit

Reactivity


Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications


Antibody used in WB on Human Muscle at 0...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ALKBH3.

Target

ALKBH3

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ALKBH3 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

ALKBH3 antibody can be used for detection of ALKBH3 by ELISA at 1:1562500. ALKBH3 antibody can be used for detection of ALKBH3 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

33 kDa

Uniprot ID
[Q96Q83](#)
NCBI
[NP_631917](#)
Dilution Range

ALKBH3 antibody can be used for detection of ALKBH3 by ELISA at 1:1562500. ALKBH3 antibody can be used for detection of ALKBH3 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.