
Product Datasheet

DHDH Antibody (orb1244445)

Description

DHDH Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody used in WB on Human Liver at 0....

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human DHDH.

Target

DHDH

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store DHDH antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

DHDH antibody can be used for detection of DHDH by ELISA at 1:1562500. DHDH antibody can be used for detection of DHDH by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

36 kDa

Uniprot ID
[Q9UQ10](#)
NCBI
[NP_055290](#)
Dilution Range

DHDH antibody can be used for detection of DHDH by ELISA at 1:1562500. DHDH antibody can be used for detection of DHDH by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.