
Product Datasheet

ARCN1 Antibody (orb1244509)

Description

ARCNI Antibody

Species/Host

Rabbit

Reactivity


Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications


Antibody used in WB on Human HepG2 at 0....

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ARCNI.

Target

ARCNI

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ARCNI antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

ARCNI antibody can be used for detection of ARCNI by ELISA at 1:312500. ARCNI antibody can be used for detection of ARCNI by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

57 kDa

Uniprot ID
[P48444](#)
NCBI
[NP_001646](#)
Dilution Range

ARCNI antibody can be used for detection of ARCNI by ELISA at 1:312500. ARCNI antibody can be used for detection of ARCNI by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.