

Product Datasheet

PAOX Antibody (orb1244569)

Description

PAOX Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PAOX.

Target

PAOX

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store PAOX antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

PAOX antibody can be used for detection of PAOX by ELISA at 1:62500. PAOX antibody can be used for detection of PAOX by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

25 kDa

Uniprot ID

[Q6QHF9](#)

NCBI

[AAS64381.1](#)

Dilution Range

PAOX antibody can be used for detection of PAOX by ELISA at 1:62500. PAOX antibody can be used for detection of PAOX by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

90kDa

40kDa

31kDa

22kDa

Antibody used in WB on Human Muscle at 0...

PAOX

90

65

40

34

22

Rabbit Anti-PAOX
Sample Type: Human Fetal Brain
Antibody Concentration: 1µg/mL

Antibody used in WB on Hum. Fetal Brain ...