
Product Datasheet

RDH12 Antibody (orb1244681)

Description	RDH12 Antibody	910 ₀
Species/Host	Rabbit	410 ₀
Reactivity	Human	310 ₀
Conjugation	Unconjugated	210 ₀
Tested Applications	ELISA, WB	Antibody used in WB on Transfected 293T ...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human RDH12.	
Target	RDH12	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store RDH12 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	RDH12 antibody can be used for detection of RDH12 by ELISA at 1:312500. RDH12 antibody can be used for detection of RDH12 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	35 kDa	
Uniprot ID	Q96NR8	
NCBI	NP_689656	
Dilution Range	RDH12 antibody can be used for detection of RDH12 by ELISA at 1:312500. RDH12 antibody can be used for detection of RDH12 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	