



Product Datasheet

SIRT1 antibody (orb12448)


Description	Rabbit polyclonal antibody to SIRT1
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic peptide of human SIRT1(NP_036370.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	120KDa
Uniprot ID	Q96EB6
Dilution Range	WB: 1:200-1:500, IHC-P: 1:100-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Immunohistochemical staining of human to...



Immunohistochemical staining of human ma...