


Product Datasheet

P2RY1 Antibody (orb1244808)


Description

P2RY1 Antibody

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human P2RY12.
Target	P2RY1
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store P2RY12 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	P2RY12 antibody can be used for detection of P2RY12 by ELISA at 1:312500. P2RY12 antibody can be used for detection of P2RY12 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	39 kDa
Uniprot ID	P47900
NCBI	NP_795345
Dilution Range	P2RY12 antibody can be used for detection of P2RY12 by ELISA at 1:312500. P2RY12 antibody can be used for detection of P2RY12 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human Fetal Muscl...



Antibody used in WB on Human 293T at 1 u...