

Product Datasheet

PGAP3 Antibody (orb1244921)

Description	PGAP3 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	Antibody used in WB on Human HeLa at 0.2...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PERLD1.	
Target	PGAP3	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store PERLD1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	PERLD1 antibody can be used for detection of PERLD1 by ELISA at 1:12500. PERLD1 antibody can be used for detection of PERLD1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	36 kDa	
Uniprot ID	Q96FM1	
NCBI	NP_219487	
Dilution Range	PERLD1 antibody can be used for detection of PERLD1 by ELISA at 1:12500. PERLD1 antibody can be used for detection of PERLD1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	