

Product Datasheet

ORAI2 Antibody (orb1244929)

Description	ORAI2 Antibody	910_
Species/Host	Rabbit	400_
Reactivity	Human	310_
Conjugation	Unconjugated	210_
Tested Applications	ELISA, WB	Antibody used in WB on Human MCF-7 at 0....
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ORAI2.	
Target	ORAI2	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store ORAI2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	ORAI2 antibody can be used for detection of ORAI2 by ELISA at 1:1562500. ORAI2 antibody can be used for detection of ORAI2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	28 kDa	
Uniprot ID	Q96SN7	
NCBI	NP_116220	
Dilution Range	ORAI2 antibody can be used for detection of ORAI2 by ELISA at 1:1562500. ORAI2 antibody can be used for detection of ORAI2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	