
Product Datasheet

SFXN3 Antibody (orb1244943)

Description	SFXN3 Antibody	90k
Species/Host	Rabbit	40k
Reactivity	Human, Mouse, Rat	20k
Conjugation	Unconjugated	20k
Tested Applications	ELISA, WB	Antibody used in WB on Human CCRF-CEM at...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SFXN3.	
Target	SFXN3	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store SFXN3 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	SFXN3 antibody can be used for detection of SFXN3 by ELISA at 1:1562500. SFXN3 antibody can be used for detection of SFXN3 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	36 kDa	
Uniprot ID	A6NCZ6	
NCBI	NP_112233	
Dilution Range	SFXN3 antibody can be used for detection of SFXN3 by ELISA at 1:1562500. SFXN3 antibody can be used for detection of SFXN3 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	