

Product Datasheet

EDAR Antibody (orb1244965)


Description	EDAR Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human EDAR.
Target	EDAR
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store EDAR antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	EDAR antibody can be used for detection of EDAR by ELISA at 1:12500. EDAR antibody can be used for detection of EDAR by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	48 kDa
Uniprot ID	Q9UNE0
NCBI	NP_071731
Dilution Range	EDAR antibody can be used for detection of EDAR by ELISA at 1:12500. EDAR antibody can be used for detection of EDAR by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human Placenta at...



Antibody used in WB on Hum. Fetal Lung a...



Antibody used in WB on Hum. Fetal Muscle...