

Product Datasheet

JPH3 Antibody (orb1244989)

Description	JPH3 Antibody	90kDa
Species/Host	Rabbit	40kDa
Reactivity	Human, Mouse, Rat	31kDa
Conjugation	Unconjugated	22kDa
Tested Applications	ELISA, WB	Antibody used in WB on Human HepG2 at 0....
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human JPH3.	
Target	JPH3	Antibody used in WB on HumanBrain at 1 u...
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store JPH3 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	JPH3 antibody can be used for detection of JPH3 by ELISA at 1:312500. JPH3 antibody can be used for detection of JPH3 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	81 kDa	
Uniprot ID	Q8WXH2	
NCBI	NP_065706	
Dilution Range	JPH3 antibody can be used for detection of JPH3 by ELISA at 1:312500. JPH3 antibody can be used for detection of JPH3 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	