

Product Datasheet

PCDHA5 Antibody (orb1245010)

Description	PCDHA5 Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	Antibody used in WB on Human HeLa at 0.2...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PCDHA5.	
Target	PCDHA5	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store PCDHA5 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	PCDHA5 antibody can be used for detection of PCDHA5 by ELISA at 1:312500. PCDHA5 antibody can be used for detection of PCDHA5 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	99 kDa	
Uniprot ID	Q9Y5H7	
NCBI	NP_061731	
Dilution Range	PCDHA5 antibody can be used for detection of PCDHA5 by ELISA at 1:312500. PCDHA5 antibody can be used for detection of PCDHA5 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	