
Product Datasheet

GALNT18 Antibody (orb1245021)

Description

GALNT18 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

 Antibody
used in WB
on Human
721_B at 0....

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human GALNTL4.

Target

GALNT18

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store GALNTL4 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

GALNTL4 antibody can be used for detection of GALNTL4 by ELISA at 1:62500. GALNTL4 antibody can be used for detection of GALNTL4 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

69 kDa

Uniprot ID
[Q6P9A2](#)
NCBI
[NP_940918](#)
Dilution Range

GALNTL4 antibody can be used for detection of GALNTL4 by ELISA at 1:62500. GALNTL4 antibody can be used for detection of GALNTL4 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.