
Product Datasheet

CAMLG Antibody (orb1245110)

Description	CAMLG Antibody	
Species/Host	Rabbit	80%
Reactivity	Human	40%
Conjugation	Unconjugated	30%
Tested Applications	ELISA, WB	20%
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CAMLG.	Antibody used in WB on Human Liver at 0....
Target	CAMLG	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store CAMLG antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	CAMLG antibody can be used for detection of CAMLG by ELISA at 1:12500. CAMLG antibody can be used for detection of CAMLG by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	33 kDa	
Uniprot ID	P49069	
NCBI	NP_001736	
Dilution Range	CAMLG antibody can be used for detection of CAMLG by ELISA at 1:12500. CAMLG antibody can be used for detection of CAMLG by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	