
Product Datasheet

MEIS3 Antibody (orb1245155)

Description	MEIS3 Antibody	
Species/Host	Rabbit	410
Reactivity	Human, Mouse, Rat	310
Conjugation	Unconjugated	210
Tested Applications	ELISA, WB	Antibody used in WB on Human Liver at 0....
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human MEIS3.	
Target	MEIS3	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store MEIS3 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	MEIS3 antibody can be used for detection of MEIS3 by ELISA at 1:62500. MEIS3 antibody can be used for detection of MEIS3 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	39 kDa	
Uniprot ID	Q99687	
NCBI	NP_001009813	
Dilution Range	MEIS3 antibody can be used for detection of MEIS3 by ELISA at 1:62500. MEIS3 antibody can be used for detection of MEIS3 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	