

Product Datasheet

KLRA1 Antibody (orb1245245)

Description	KLRA1 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	 <p data-bbox="1106 497 1230 622">Antibody used in WB on Human Liver at 0....</p>
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human KLRA1.	
Target	KLRA1	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store KLRA1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	KLRA1 antibody can be used for detection of KLRA1 by ELISA at 1:312500. KLRA1 antibody can be used for detection of KLRA1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	25 kDa	
Uniprot ID	Q9UNGO	
NCBI	NP_006602	
Dilution Range	KLRA1 antibody can be used for detection of KLRA1 by ELISA at 1:312500. KLRA1 antibody can be used for detection of KLRA1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	