
Product Datasheet

ST3GAL4 Antibody (orb1245252)

Description

ST3GAL4 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications
Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ST3GAL4.

Target

ST3GAL4

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ST3GAL4 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

ST3GAL4 antibody can be used for detection of ST3GAL4 by ELISA at 1:62500. ST3GAL4 antibody can be used for detection of ST3GAL4 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

37 kDa

Uniprot ID
[Q6IBE6](#)
NCBI
[NP_006269](#)
Dilution Range


ST3GAL4 antibody can be used for detection of ST3GAL4 by ELISA at 1:62500. ST3GAL4 antibody can be used for detection of ST3GAL4 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date


12 months from date of receipt.



Antibody used in WB on Human 721_B at 0...



Antibody used in WB on Human HeLa at 1 u...



Antibody used in WB on Human MCF7 at 1 u...