

Product Datasheet

SGCE Antibody (orb1245273)


Description	SGCE Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SGCE.
Target	SGCE
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store SGCE antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	SGCE antibody can be used for detection of SGCE by ELISA at 1:312500. SGCE antibody can be used for detection of SGCE by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	52 kDa
Uniprot ID	B5MDA7
NCBI	NP_001092871
Dilution Range	SGCE antibody can be used for detection of SGCE by ELISA at 1:312500. SGCE antibody can be used for detection of SGCE by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in IHC on human skeletal m...



Antibody used in WB on Human 721_B at 0...



Antibody used in IHC on Mouse Cerebellum...