
Product Datasheet

OXSM Antibody (orb1245286)

Description	OXSM Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human OXSM.	Antibody
Target	OXSM	used in WB
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	on
Form/Appearance	Liquid	Transfected
Concentration	batch dependent	293T ...
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store OXSM antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	OXSM antibody can be used for detection of OXSM by ELISA at 1:312500. OXSM antibody can be used for detection of OXSM by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	49 kDa	
Uniprot ID	Q9NWU1	
NCBI	NP_060367	
Dilution Range	OXSM antibody can be used for detection of OXSM by ELISA at 1:312500. OXSM antibody can be used for detection of OXSM by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	