
Product Datasheet

PSMG1 Antibody (orb1245328)

Description	PSMG1 Antibody	
Species/Host	Rabbit	80%
Reactivity	Human	40%
Conjugation	Unconjugated	20%
Tested Applications	ELISA, WB	20%
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PSMG1.	Antibody used in WB on Human Heart at 0....
Target	PSMG1	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store PSMG1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	PSMG1 antibody can be used for detection of PSMG1 by ELISA at 1:62500. PSMG1 antibody can be used for detection of PSMG1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	33 kDa	
Uniprot ID	O95456	
NCBI	NP_003711	
Dilution Range	PSMG1 antibody can be used for detection of PSMG1 by ELISA at 1:62500. PSMG1 antibody can be used for detection of PSMG1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	