
Product Datasheet

PRDM15 Antibody (orb1245329)

Description

PRDM15 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

 Antibody
used in WB
on Human
293T at 0.2...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PRDM15.

Target

PRDM15

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store PRDM15 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

PRDM15 antibody can be used for detection of PRDM15 by ELISA at 1:312500. PRDM15 antibody can be used for detection of PRDM15 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

134 kDa

Uniprot ID
[P57071](#)
NCBI
[NP_001035514](#)
Dilution Range

PRDM15 antibody can be used for detection of PRDM15 by ELISA at 1:312500. PRDM15 antibody can be used for detection of PRDM15 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.