
Product Datasheet

SLCO2B1 Antibody (orb1245496)

Description

SLCO2B1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SLCO2B1.

Immunogen
Target

SLCO2B1

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store SLCO2B1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

SLCO2B1 antibody can be used for detection of SLCO2B1 by ELISA at 1:62500. SLCO2B1 antibody can be used for detection of SLCO2B1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

77 kDa

Uniprot ID
[O94956](https://www.uniprot.org/uniprot/Q94956)
NCBI
[NP_009187](https://www.ncbi.nlm.nih.gov/nuccore/NP_009187)
Dilution Range

SLCO2B1 antibody can be used for detection of SLCO2B1 by ELISA at 1:62500. SLCO2B1 antibody can be used for detection of SLCO2B1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



40kDa

30kDa

20kDa

Antibody used in WB on Human HepG2 cells...



90

65

40

25

22

Rabbit Anti-SLCO2B1

Sample Type Human Fetal Lung

Antibody Concentration 1µg/ml

Antibody used in WB on Hum. Fetal Lung a...



90

65

40

25

22

Rabbit Anti-SLCO2B1

Sample Type Human Adult Placenta

Antibody Concentration 1µg/ml

Antibody used in WB on Hum. Adult Placen...