

## Product Datasheet

### SLC10A7 Antibody (orb1245511)

**Description**

SLC10A7 Antibody

**Species/Host**

Rabbit

**Reactivity**


Human, Mouse, Rat

**Conjugation**

Unconjugated

**Tested**

ELISA, WB

**Applications**


Antibody used in WB on Human HepG2 at 0....

**Immunogen**

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SLC10A7.

**Target**

SLC10A7

**Preservatives**

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

**Form/Appearance**

Liquid

**Concentration**

batch dependent

**Storage**

For short periods of storage (days) store at 4°C. For longer periods of storage, store SLC10A7 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

**Note**

For research use only

**Application notes**

SLC10A7 antibody can be used for detection of SLC10A7 by ELISA at 1:1562500. SLC10A7 antibody can be used for detection of SLC10A7 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

**Clonality**

Polyclonal

**MW**

37 kDa

**Uniprot ID**
[Q0GE19](#)
**NCBI**
[NP\\_001025169](#)
**Dilution Range**

SLC10A7 antibody can be used for detection of SLC10A7 by ELISA at 1:1562500. SLC10A7 antibody can be used for detection of SLC10A7 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

**Expiration Date**

12 months from date of receipt.