

Product Datasheet

HAAO Antibody (orb1245643)

Description	HAAO Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	Antibody used in WB on Human 293T at 0.2...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human HAAO.	
Target	HAAO	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store HAAO antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	HAAO antibody can be used for detection of HAAO by ELISA at 1:312500. HAAO antibody can be used for detection of HAAO by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	32 kDa	
Uniprot ID	P46952	
NCBI	NP_036337	
Dilution Range	HAAO antibody can be used for detection of HAAO by ELISA at 1:312500. HAAO antibody can be used for detection of HAAO by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	