


Product Datasheet

GPR6 Antibody (orb1245655)

Description	GPR6 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human GPR6.	Antibody used in WB on Human Placenta at...
Target	GPR6	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	Antibody used in IHC on Human Placenta a...
Form/Appearance	Liquid	
Concentration	batch dependent	Antibody used in IHC on Human Placenta a...
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store GPR6 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	Antibody used in IHC on Human Placenta a...
Note	For research use only	Antibody used in IHC on Human Placenta a...
Application notes	GPR6 antibody can be used for detection of GPR6 by ELISA at 1:312500. GPR6 antibody can be used for detection of GPR6 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	Antibody used in IHC on Human Placenta a...
Clonality	Polyclonal	Antibody used in IHC on Human Placenta a...
MW	38 kDa	Antibody used in IHC on Human Placenta a...
Uniprot ID	P46095	Antibody used in IHC on Human Placenta a...
NCBI	NP_005275	Antibody used in IHC on Human Placenta a...
Dilution Range	GPR6 antibody can be used for detection of GPR6 by ELISA at 1:312500. GPR6 antibody can be used for detection of GPR6 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	Antibody used in IHC on Human Placenta a...
Expiration Date	12 months from date of receipt.	Antibody used in IHC on Human Placenta a...