

Product Datasheet

UPP1 Antibody (orb1245663)

Description		UPP1 Antibody	
Species/Host	Rabbit		
Reactivity	Human		
Conjugation	Unconjugated		
Tested Applications	ELISA, WB		Antibody used in WB on Human 293T at 0.2...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human UPP1.		
Target	UPP1		
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.		
Form/Appearance	Liquid		
Concentration	batch dependent		
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store UPP1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.		
Note	For research use only		
Application notes	UPP1 antibody can be used for detection of UPP1 by ELISA at 1:62500. UPP1 antibody can be used for detection of UPP1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.		
Clonality	Polyclonal		
MW	34 kDa		
Uniprot ID	Q16831		
NCBI	NP_003355		
Dilution Range	UPP1 antibody can be used for detection of UPP1 by ELISA at 1:62500. UPP1 antibody can be used for detection of UPP1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.		
Expiration Date	12 months from date of receipt.		