
Product Datasheet

LAIR2 Antibody (orb1245666)

Description	LAIR2 Antibody	
Species/Host	Rabbit	80kD
Reactivity	Human	41kD
Conjugation	Unconjugated	31kD
Tested Applications	ELISA, WB	21kD
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human LAIR2.	11kD
Target	LAIR2	Antibody used in WB on Human MCF-7 at 0....
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store LAIR2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	LAIR2 antibody can be used for detection of LAIR2 by ELISA at 1:312500. LAIR2 antibody can be used for detection of LAIR2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	16 kDa	
Uniprot ID	Q6ISS4	
NCBI	NP_002279	
Dilution Range	LAIR2 antibody can be used for detection of LAIR2 by ELISA at 1:312500. LAIR2 antibody can be used for detection of LAIR2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	