

Product Datasheet

Caspase 3 antibody (orb12457)

Description Rabbit polyclonal antibody to Caspase 3

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ICC, IF, IHC, WB

Immunogen A synthetic peptide of human Caspase-3(NP_116786.1).

Preservatives PBS with 0.05% proclin300,50% glycerol,pH7.3.

Storage Store at -20°C. Avoid freeze / thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 17KDa/19KDa/32KDa


Uniprot ID [P42574](#)

Dilution Range WB: 1:500-1:2000, IF/ICC: 1:50-1:200


Expiration Date 12 months from date of receipt.



Immunohistochemical staining of rat kidn...



Immunohistochemical staining of rat test...



Immunohistochemical staining of rat hear...