

---

## Product Datasheet

### SNRPD2 Antibody (orb1245778)

**Description**

SNRPD2 Antibody

**Species/Host**

Rabbit

**Reactivity**

Human, Mouse, Rat

**Conjugation**

Unconjugated

**Tested**

ELISA, WB

**Applications**

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SNRPD2.

**Immunogen**
**Target**

SNRPD2

**Preservatives**

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

**Form/Appearance**

Liquid

**Concentration**

batch dependent

**Storage**

For short periods of storage (days) store at 4°C. For longer periods of storage, store SNRPD2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

**Note**

For research use only

**Application notes**

SNRPD2 antibody can be used for detection of SNRPD2 by ELISA at 1:312500. SNRPD2 antibody can be used for detection of SNRPD2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

**Clonality**

Polyclonal

**MW**

13 kDa

**Uniprot ID**
[P62316](#)
**NCBI**
[NP\\_004588](#)
**Dilution Range**

SNRPD2 antibody can be used for detection of SNRPD2 by ELISA at 1:312500. SNRPD2 antibody can be used for detection of SNRPD2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

**Expiration Date**

12 months from date of receipt.







Antibody used in WB on Human Jurkat 0.2-...










Antibody used in WB on Human MCF7 at 1 u...









Antibody used in WB on Human 293T at 1 u...