

Product Datasheet

TUBA4A Antibody (orb1245810)

Description

TUBA4A Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human TUBA4A.

Target

TUBA4A

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store TUBA4A antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

TUBA4A antibody can be used for detection of TUBA4A by ELISA at 1:12500. TUBA4A antibody can be used for detection of TUBA4A by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

50 kDa

Uniprot ID

[P68366](#)

NCBI

[NP_005991](#)


Dilution Range

TUBA4A antibody can be used for detection of TUBA4A by ELISA at 1:12500. TUBA4A antibody can be used for detection of TUBA4A by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

TUBA4A




Antibody used in WB on Mouse Brain at 2 ...



Antibody used in WB on Human Placenta at...

TUBA4A



Antibody-RED
DAPI-Blue

(See IHC Data and Customer Feedback for more information)

Antibody used in IHC on Mouse brain stem...