
Product Datasheet

GRHL2 Antibody (orb1245864)

Description

GRHL2 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications
Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human GRHL2.

Target

GRHL2

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store GRHL2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

GRHL2 antibody can be used for detection of GRHL2 by ELISA at 1:12500. GRHL2 antibody can be used for detection of GRHL2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW


71 kDa

Uniprot ID
[Q6ISB3](#)
NCBI
[NP_079191](#)
Dilution Range


GRHL2 antibody can be used for detection of GRHL2 by ELISA at 1:12500. GRHL2 antibody can be used for detection of GRHL2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC.



Antibody used in IHC.



Antibody used in WB on Human 721_B at 1 ...