

Product Datasheet

ZBTB3 Antibody (orb1245866)

Description	ZBTB3 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	Antibody used in WB on Human HeLa at 0.2...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ZBTB3.	
Target	ZBTB3	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store ZBTB3 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	Antibody used in WB on Hum. Adult Placen...
Note	For research use only	
Application notes	ZBTB3 antibody can be used for detection of ZBTB3 by ELISA at 1:62500. ZBTB3 antibody can be used for detection of ZBTB3 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	62 kDa	
Uniprot ID	Q9H5J0	
NCBI	NP_079060	
Dilution Range	ZBTB3 antibody can be used for detection of ZBTB3 by ELISA at 1:62500. ZBTB3 antibody can be used for detection of ZBTB3 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	



Antibody used in WB on Human HeLa at 0.2...



Antibody used in WB on Hum. Adult Placen...