

Product Datasheet

HOXC8 Antibody (orb1245871)

Description

HOXC8 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human HOXC8.

Target

HOXC8

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store HOXC8 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

HOXC8 antibody can be used for detection of HOXC8 by ELISA at 1:312500. HOXC8 antibody can be used for detection of HOXC8 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

28 kDa

Uniprot ID

[P31273](#)

NCBI


[NP_073149](#)

Dilution Range


HOXC8 antibody can be used for detection of HOXC8 by ELISA at 1:312500. HOXC8 antibody can be used for detection of HOXC8 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date


12 months from date of receipt.



Antibody used in IHC on Human small inte...



Antibody used in IHC on Human Colon at 5...



Antibody used in WB on Human Fetal Muscl...