

---

## Product Datasheet

### KLHL1 Antibody (orb1245887)

Description	KLHL1 Antibody	910k <sub>Da</sub>
<b>Species/Host</b>	Rabbit	410k <sub>Da</sub>
<b>Reactivity</b>	Human, Mouse	310k <sub>Da</sub>
<b>Conjugation</b>	Unconjugated	210k <sub>Da</sub>
<b>Tested Applications</b>	ELISA, WB	Antibody used in WB on Human 293T at 0.2...
<b>Immunogen</b>	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human KLHL1.	KLHL1
<b>Target</b>	KLHL1	90
<b>Preservatives</b>	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	65
<b>Form/Appearance</b>	Liquid	40
<b>Concentration</b>	batch dependent	31
<b>Storage</b>	For short periods of storage (days) store at 4°C. For longer periods of storage, store KLHL1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	22
<b>Note</b>	For research use only	Antibody used in WB on Human 721_B at 1 ...
<b>Application notes</b>	KLHL1 antibody can be used for detection of KLHL1 by ELISA at 1:62500. KLHL1 antibody can be used for detection of KLHL1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	KLHL1
<b>Clonality</b>	Polyclonal	90
<b>MW</b>	83 kDa	65
<b>Uniprot ID</b>	<a href="#">Q9NR64</a>	40
<b>NCBI</b>	<a href="#">NP_065917</a>	31
<b>Dilution Range</b>	KLHL1 antibody can be used for detection of KLHL1 by ELISA at 1:62500. KLHL1 antibody can be used for detection of KLHL1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	22
<b>Expiration Date</b>	12 months from date of receipt.	