
Product Datasheet

MNT Antibody (orb1245893)

Description	MNT Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	Antibody used in WB on Human HeLa at 0.2...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human MNT.	
Target	MNT	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store MNT antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	Antibody used in WB on Hum. Fetal Lung a...
Note	For research use only	
Application notes	MNT antibody can be used for detection of MNT by ELISA at 1:1562500. MNT antibody can be used for detection of MNT by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	62 kDa	
Uniprot ID	Q99583	Antibody used in WB on Human HepG2 at 1 ...
NCBI	NP_064706	
Dilution Range	MNT antibody can be used for detection of MNT by ELISA at 1:1562500. MNT antibody can be used for detection of MNT by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	

