
Product Datasheet

GRIK5 Antibody (orb1246022)

Description

GRIK5 Antibody

Species/Host

Rabbit

Reactivity


Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications


Antibody used in WB on Human Jurkat 0.2-...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human GRIK5.

Target

GRIK5

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store GRIK5 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

GRIK5 antibody can be used for detection of GRIK5 by ELISA at 1:62500. GRIK5 antibody can be used for detection of GRIK5 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

108 kDa

Uniprot ID
[Q16478](#)
NCBI
[NP_002079](#)
Dilution Range

GRIK5 antibody can be used for detection of GRIK5 by ELISA at 1:62500. GRIK5 antibody can be used for detection of GRIK5 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.