

Product Datasheet

CHRFAM7A Antibody (orb1246143)

Description

CHRFAM7A Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CHRFAM7A.

Target

CHRFAM7A

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store CHRFAM7A antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

CHRFAM7A antibody can be used for detection of CHRFAM7A by ELISA at 1:62500. CHRFAM7A antibody can be used for detection of CHRFAM7A by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

35 kDa

Uniprot ID

[Q8IUZ4](#)

NCBI

[NP_683709](#)


Dilution Range

CHRFAM7A antibody can be used for detection of CHRFAM7A by ELISA at 1:62500. CHRFAM7A antibody can be used for detection of CHRFAM7A by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.


Expiration Date

12 months from date of receipt.

Antibody used in WB on Human 293T at 0.2...



Antibody used in WB on Hum. Fetal Muscle...



Antibody used in WB on Human 721_B at 1 ...