

Product Datasheet

VDAC3 Antibody (orb1246156)


Description	VDAC3 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human VDAC3.
Target	VDAC3
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store VDAC3 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	VDAC3 antibody can be used for detection of VDAC3 by ELISA at 1:62500. VDAC3 antibody can be used for detection of VDAC3 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	31 kDa
Uniprot ID	Q9Y277
NCBI	NP_005653
Dilution Range	VDAC3 antibody can be used for detection of VDAC3 by ELISA at 1:62500. VDAC3 antibody can be used for detection of VDAC3 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human Heart at 0....



Antibody used in WB on Hum. Fetal Muscle...



Antibody used in WB on Hum. Fetal Heart ...