
Product Datasheet

RBBP5 Antibody (orb1246213)

Description	RBBP5 Antibody	RBBP5
Species/Host	Rabbit	40kDa
Reactivity	Human, Mouse, Rat	31kDa
Conjugation	Unconjugated	22kDa
Tested Applications	ELISA, WB	Antibody used in WB on Human HepG2 at 0....
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human RBBP5.	RBBP5
Target	RBBP5	90
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	65
Form/Appearance	Liquid	40
Concentration	batch dependent	31
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store RBBP5 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	22
Note	For research use only	Rabbit Anti-RBBP5
Application notes	RBBP5 antibody can be used for detection of RBBP5 by ELISA at 1:12500. RBBP5 antibody can be used for detection of RBBP5 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	Sample Type: Human Adult Placenta
Clonality	Polyclonal	Antibody Concentration: 1mg/ml
MW	59 kDa	RBBP5
Uniprot ID	Q15291	90
NCBI	NP_005048	65
Dilution Range	RBBP5 antibody can be used for detection of RBBP5 by ELISA at 1:12500. RBBP5 antibody can be used for detection of RBBP5 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	40
Expiration Date	12 months from date of receipt.	31
		22
		Rabbit Anti-RBBP5
		Sample Type: Human MCF7
		Antibody Concentration: 1mg/ml
		Antibody used in WB on Human MCF7 at 1 u...