

Product Datasheet

RXRA Antibody (orb1246222)


Description

RXRA Antibody


Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human RXRA.
Target	RXRA
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store RXRA antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	RXRA antibody can be used for detection of RXRA by ELISA at 1:312500. RXRA antibody can be used for detection of RXRA by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	51 kDa
Uniprot ID	P19793
NCBI	NP_002948
Dilution Range	RXRA antibody can be used for detection of RXRA by ELISA at 1:312500. RXRA antibody can be used for detection of RXRA by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody
used in WB
on
Transfected
293T ...



Antibody
used in WB
on Hum.
Fetal Heart ...



Antibody
used in WB
on Hum.
Fetal Lung
a...