
Product Datasheet

RGS6 Antibody (orb1246223)

Description	RGS6 Antibody	EVID_1
Species/Host	Rabbit	4100_
Reactivity	Human, Rat	3100_
Conjugation	Unconjugated	2100_
Tested Applications	ELISA, WB	Antibody used in WB on Human Jurkat 0.2-...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human RGS6.	
Target	RGS6	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	RGS6
Form/Appearance	Liquid	90_
Concentration	batch dependent	65_
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store RGS6 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	40_
Note	For research use only	28_
Application notes	RGS6 antibody can be used for detection of RGS6 by ELISA at 1:12500. RGS6 antibody can be used for detection of RGS6 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	22_
Clonality	Polyclonal	Rabbit Anti-RGS6 Sample Type: Human T21_B Antibody Concentration: 1µg/mL
MW	54 kDa	
Uniprot ID	P49758	
NCBI	NP_004287	
Dilution Range	RGS6 antibody can be used for detection of RGS6 by ELISA at 1:12500. RGS6 antibody can be used for detection of RGS6 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	