
Product Datasheet

TAF2 Antibody (orb1246266)

Description	TAF2 Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human TAF2.	Antibody used in WB on Human 293T at 0.2...
Target	TAF2	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store TAF2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	TAF2 antibody can be used for detection of TAF2 by ELISA at 1:12500. TAF2 antibody can be used for detection of TAF2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	137 kDa	
Uniprot ID	Q6P1X5	
NCBI	NP_003175	
Dilution Range	TAF2 antibody can be used for detection of TAF2 by ELISA at 1:12500. TAF2 antibody can be used for detection of TAF2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	