
Product Datasheet

JMJD4 Antibody (orb1246273)

Description	JMJD4 Antibody	EVID_
Species/Host	Rabbit	4100_
Reactivity	Human	3100_
Conjugation	Unconjugated	2100_
Tested Applications	ELISA, WB	Antibody used in WB on Human Heart at 0...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human JMJD4.	JMJD4
Target	JMJD4	90_
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	65_
Form/Appearance	Liquid	40_
Concentration	batch dependent	29_
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store JMJD4 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	22_
Note	For research use only	Antibody used in WB on Human 721_B at 1 ...
Application notes	JMJD4 antibody can be used for detection of JMJD4 by ELISA at 1:12500. JMJD4 antibody can be used for detection of JMJD4 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	JMJD4
Clonality	Polyclonal	90_
MW	52 kDa	65_
Uniprot ID	Q9H9V9	40_
NCBI	NP_075383	29_
Dilution Range	JMJD4 antibody can be used for detection of JMJD4 by ELISA at 1:12500. JMJD4 antibody can be used for detection of JMJD4 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	22_
Expiration Date	12 months from date of receipt.	Antibody used in WB on Hum. Fetal Liver ...