

Product Datasheet

NKX3-2 Antibody (orb1246279)

Description

NKX3-2 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody used in WB on Human Jurkat 0.2-...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human NKX3-2.

Target

NKX3-2

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store NKX3-2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Antibody used in WB on Hum. Fetal Muscle...

Note

For research use only

Application notes

NKX3-2 antibody can be used for detection of NKX3-2 by ELISA at 1:62500. NKX3-2 antibody can be used for detection of NKX3-2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

35 kDa

Uniprot ID
[P78367](#)
NCBI
[NP_001180](#)
Dilution Range

NKX3-2 antibody can be used for detection of NKX3-2 by ELISA at 1:62500. NKX3-2 antibody can be used for detection of NKX3-2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

