
Product Datasheet

LHX4 Antibody (orb1246336)

Description	LHX4 Antibody	Bio
Species/Host	Rabbit	Bio
Reactivity	Human, Mouse, Rat	Bio
Conjugation	Unconjugated	Bio
Tested Applications	ELISA, WB	Bio
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human LHX4.	Antibody used in WB on Human 721_B at 0....
Target	LHX4	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store LHX4 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	LHX4 antibody can be used for detection of LHX4 by ELISA at 1:12500. LHX4 antibody can be used for detection of LHX4 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	43 kDa	
Uniprot ID	P53776	
NCBI	NP_203129	
Dilution Range	LHX4 antibody can be used for detection of LHX4 by ELISA at 1:12500. LHX4 antibody can be used for detection of LHX4 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	