




Product Datasheet

GSC Antibody (orb1246354)

Description	GSC Antibody	
Species/Host	Rabbit	 
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human GSC.	Antibody used in IHC on Human Tonsil at ...
Target	GSC	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store GSC antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	Antibody used in IHC on Human Testis at ...
Note	For research use only	
Application notes	GSC antibody can be used for detection of GSC by ELISA at 1:12500. GSC antibody can be used for detection of GSC by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	28 kDa	
Uniprot ID	P56915	Antibody used in WB on Human 721_B at 0....
NCBI	NP_776248	
Dilution Range	GSC antibody can be used for detection of GSC by ELISA at 1:12500. GSC antibody can be used for detection of GSC by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	